







CERTIFICATE OF CONFORMITY

EU - Medical devices Regulation (MDR)-
No. CIC200309003-ZJ-ZS

Applicant : BeijingYimed TechnologyCo. ,Ltd.
Address : Room 504,5th floor,Nailunbuilding,No.1,LiuliBridge,Fengtai District,Beijing
Manufacturer : BeijingYimedTechnologyCo. ,Ltd.
Address : Room 504,5th floor,Nailunbuilding,No.1,LiuliBridge,Fengtai District,Beijing
Product : Surgical Mask
Trade Name : N/A
Model(s) : YMD-188
Test Report No. : CIC200309003-ZJ
Test Standards : EN 14683:2019

The tests that base on the above designated product Complies with the essential requirements of Regulation 2017/745/EU.

The test results apply only to the particular sample tested and to the specific tests carried out.

Specifically to the sample investigated in our test reference number only.

The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation.

Other relevant Directives have to be observed..



Hale Zhao (Chief engineer)
March 11, 2020

Shenzhen CIC Testing Technology Co., Ltd

Add: No.102, Floor 1, Building A, No. 26 of Xinhe Road, Xinqiao Community, Xinqiao Street,
Baoan District, Shenzhen, Guangdong, China.

Web: [http:// www.soncap-saso.com](http://www.soncap-saso.com) Tel: 0755-23593247 E-mail: admin@cic-test.com



TEST REPORT

MEASUREMENT AND TEST REPORT For

BeijingYimed TechnologyCo. ,Ltd.

Room 504,5th floor,Nailunbuilding,No.1,LiuliBridge,FengtaiDistrict,Beijing

Models: YMD-188

March 11, 2020

This Report Concerns: <input checked="" type="checkbox"/> Original Report	Equipment Type: Surgical Mask
Test Engineer:	<i>Dana Huang</i>
Report Number:	CIC200309003-ZJ
Test Date:	March 09-10, 2020
Reviewed By:	 <i>Hong Zhang</i>
Prepared By:	Shenzhen CIC Testing Technology Co., Ltd. No.102, Floor 1, Building A, No. 26 of Xinhe Road, Xinqiao Community, Xinqiao Street, Baoan District, Shenzhen, Guangdong, China. Tel: 0755-23593247 Fax: +86-0755-23208027 Web: http:// www.soncap-saso.com

Note: This test report is limited to the above client company and the product model only. It may not be duplicated without prior written consent of Shenzhen CIC Testing Technology Co., Ltd.

TEST REPORT EN 14683 Medical face masks - Requirements and test methods	
Report reference No. :	CIC200309003-ZJ
Tested by (+signature)..... :	Dana Huangsee cover page.....
Approved by (+ signature).. :	Hale Zhaosee cover page.....
Date of issue..... :	March 11, 2020
Testing laboratory	
Name..... :	Shenzhen CIC Testing Technology Co., Ltd
Address..... :	No.102, Floor 1, Building A, No. 26 of Xinhe Road, Xinqiao Community, Xinqiao Street, Baoan District, Shenzhen, Guangdong, China.
Test location	(Same as above)
Client	
Name..... :	BeijingYimed TechnologyCo. ,Ltd.
Address..... :	Room 504,5th floor,Nailunbuilding,No.1,LiuliBridge,FengtaiDistrict,Beijing
Test specification	
Standard..... :	EN 14683:2019
Test procedure..... :	CE- MDR
Procedure deviation..... :	N.A.
Non-standard test method.. :	N.A.
Test Report Form No.: EN_14683	
Test item	
Description	Surgical Mask
Model No. :	YMD-188
Trade Mark..... :	N/A
Manufacturer..... :	BeijingYimed TechnologyCo. ,Ltd.
Address..... :	Room 504,5th floor,Nailunbuilding,No.1,LiuliBridge,FengtaiDistrict,Beijing
Classification..... :	Type II

Test case verdicts

Test case does not apply to the test object : N(.A.)

Test item does meet the requirement : P(ass)

Test item does not meet the requirement : F(ail)

Testing

Date of receipt of test item: March 09, 2020

Date(s) of performance of test: March 09-10, 2020

General remarks

This test report shall not be reproduced except in full without the written approval of the testing laboratory.

The test results presented in this report relate only to the item tested.

"(see remark #)" refers to a remark appended to the report.

"(see appended table)" refers to a table appended to the report.

Throughout this report a comma is used as the decimal separator.

Remark :

The EUT complies with the standard EN 14683 requirement.

Copy of marking plate:

Table 1 — Performance requirements for medical face masks


Table 1	The performance test			Result
Test	Type I a	Type II	Type IIR	Type II
Bacterial filtration efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98	99
Differential pressure (Pa/cm ²)	< 40	< 40	< 60	37.9
Splash resistance pressure (kPa)	Not required	Not required	$\geq 16,0$	--
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30	28.54
Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.				

Appendix
Photo documentation

Photo 1

View:

- ☒ front
☐ rear
☐ right side
☐ left side
☐ top
☐ bottom
☐ internal



***** End of this report *****